



# AO Infinity Foot Detox



The AO Infinity Foot Detox uses Radio Frequency to raise the body's internal temperature. This device circulates Radio Frequency vibration energy through the soles of the feet, moving up to the rest of the body. When the heat accelerates your metabolism, toxins, and waste are broken down and sent to the liver for elimination.

**Loran Swensen, Solex President** explains how the Food Detox system was developed, and provides an overview of its functions:

**VIDEO** [5 minutes]

[https://youtu.be/Esq\\_rdDkXsE?t=543](https://youtu.be/Esq_rdDkXsE?t=543)

## Benefits of Heat:

Blood circulation slows when the body is cold, and the cell metabolism weakens. Every time the body temperature drops one degree, the metabolism is reduced by about 12%, and immunity is reduced by more than 30%.

When an individual's temperature is at 95.9° Fahrenheit, cells multiply. Other viruses and bacteria may invade, and various functional disorders may occur.

On the other hand, when an individual's temperature is at 103.2° Fahrenheit, cells begin to die, and many viruses and bacteria cannot survive.

## Recommendation:

- Use once a day, twice at most, but the interval between should be more than three hours.
- The initial use should be 30 minutes at maximum.
- Each consecutive use can be 40-45 minutes long. The intensity should be adjusted to the user's comfort level every time.
- For personal use in a treatment cycle or a clinic, use once a day for ten days. Adjust settings every other session according to the user's comfort. Chronic diseases can be treated with a 30-day cycle.

## How it works:

When you place your bare feet on the plates of the device, it will generate a frequency of one megahertz. The device pulses the frequency into the body in 10 millisecond bursts. This delivery mechanism causes the body to go into a low-grade fever, raising the body temperature by one and a half degrees.

\* These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.

Jocelyne Colombe Holistic Nurse Practitioner, Quantum Living Advocate (860) 830-1115  
To purchase, <https://linkbuilder.ziplingo.com/XFUM>

### How to use the AO Infinity Foot Detox:

- Connect the power cord firmly into the power port of the AO Infinity Foot Detox device and the other end to a 100-200 volt outlet. Switch the power to 'on' on the back of the device. (If the power cord is stiff, ensure it is fully inserted).
- Position the unit so that the electrode plates are at a comfortable angle for the feet. You may use the stand to position the device at your preferred incline level.
- Place bare feet firmly on the electrode plates.
- Adjust the desired time of the screen display located on the device. The default is 30 minutes.
- Adjust the intensity of the energy output from low to more intense based on the user's comfort.
- Press the start button. (Troubleshooting if not powering on: check the fuse).

Energy will initially be concentrated at the ankle and slowly move up the legs and eventually the entire body; this may take several sessions.

Drinking plenty of room-temperature purified water during the session will help raise the body's temperature, increase blood flow, and support the release of toxins. As the energy oscillating waves penetrate deep into the body from the soles of the feet, the body temperature will increase.

#### Foot Detox Machine:

- ✓ Lightweight machine, easy to move and carry.
- ✓ Screen display where you will see the settings you select from the settings buttons below it.
- ✓ The blue lid can be lifted easily for your convenience.
- ✓ **The electrode plates can accommodate all foot sizes and stay at room temperature.**

#### Included:

**Stand:** The stand has a quick-adjust locking system that lets you adjust the angle of the machine. You can set up the stand angle that fits your personal needs. The stand also includes a supporting ridge to secure the device in place.

#### Optional and sold separately:

Use with the belt around the waist for greater results.



\* These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.

\*\* The buyer is responsible for making any adjustments necessary to their power supply to make this device function properly. Solex is not responsible for any replacements due to incorrect power supply issues.

**For purchasing**, click on this link: <https://linkbuilder.ziplingo.com/XFUM>



# AO Infinity Waistband



## What is it?

Optimize the benefits of radio frequencies from the AO Infinity Foot Detox with the AO Infinity Waistband! The band allows you to use the same one megahertz frequency, infrared heat, and vibration from the Foot Detox device to stimulate the liver, pancreas, gallbladder, and spleen.

**Dimensions:** 51 Inches x 10.5 Inches.

## Using with the Foot Bath:

The AO Infinity Foot Detox uses Radio Frequency to raise the body's internal temperature. This device circulates Radio Frequency vibration energy through the soles of the feet, moving up to the rest of the body. When the heat accelerates your metabolism, toxins, and waste are broken down and sent to the liver for elimination. This device is sold separately from the band.

## How to use:

- The waistband plugs into the back of the Foot Detox machine.
- Place it around the waist and secure it in place.
- Adjust the settings on the screen display located on the device.
- The model, intensity+, intensity-, and temperature settings are designated for the band. (These settings do not adjust the electrode plates' settings but the waistband's settings).
- The time, energy+, energy-, and start/stop settings are designated specifically for the electrode plates.
- The band is an optional accessory that can optimize the benefits of radio frequencies from the device.

## Purchase link:

<https://linkbuilder.ziplingo.com/XFUW>



## DISCOUNT

when the AO Infinity Waistband and Foot Bath are purchased together or as a Preferred Customer.



INDEPENDENT  
QUANTUM LIVING ADVOCATE

**SOLEX**

JOCELYNE COLOMBE APRN

(860)830-1115 JCDOVE.COM

AOSCAN4YOU@YAHOO.COM



\* These statements have not been evaluated by the Food and Drug Administration. The AO Infinity Foot Detox and Waistband devices are NOT medical devices and do not diagnose, prevent, treat or cure any disease.

\*\* The buyer is responsible for making any adjustments necessary to their power supply to make this device function properly. Solex is not responsible for any replacements due to incorrect power supply issues.